

Airway Management Training Through Updating Resusci Anne



Resusci Anne Airway Trainer Update Kit

With complete anatomy to the vocal cords, this airway trainer allows education in basic and intermediate airway techniques. Highly elastic and durable skin and tissue simulation provides lifelike feel and realistic movement

- **Educationally effective:** learn core skills in basic and intermediate airway management
- **Cost efficient:** update your existing Resusci Anne manikin
- **Clinically realistic:** anatomical accuracy and realism



Laerdal
helping save lives

Resusci Anne Airway Trainer Update Kits

Ordering Information:

316000: Resusci Anne Airway Trainer - Basic

316200: Resusci Anne Airway Trainer - SkillReporter

Includes:

*Airway Head, lung, stomach, Resusci Anne Simulator Eye Set,
I can airway lubricant, Directions for Use*

Technical Specifications

316000: Resusci Anne Airway Trainer - Basic

Fits the following Resusci Anne Models:

- Resusci Anne Basic: 310005, 310006, 310035
- Resusci Anne SkillGuide: 310015, 310016, 310045, 320060, 320080
- All discontinued Resusci Anne models that have disposable airways (except Resusci Anne Skillmeter).

316200: Resusci Anne Airway Trainer - SkillReporter

Fits the following Resusci Anne Models:

- Resusci Anne SkillReporter: 310025, 310055, 320070, 320090, 320025, 320055
- Resusci Anne SkillMeter: 300010

Training Techniques

- Head tilt / chin lift and jaw thrust maneuvers
- Cricoid Pressure
- Pupil Assessment (normal, constricted and dilated)
- Oropharyngeal and Nasopharyngeal Tube Insertion
- Bag-Valve-Mask Ventilation
- Laryngeal Mask Airway
- Combitube®
- Laryngeal Tube Airway



316000: Resusci Anne Airway Trainer - Basic



316200: Resusci Anne Airway Trainer - SkillReporter

For more product information: www.laerdal.com



Laerdal
helping save lives